







Product:	GentleSafe® EXTRA	
Trademark:	SAFE®	
REF:	879	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT III according to Regulation (EU) 2016/425 for personal protective equipment (PPE). Class I according to Regulation (EU) 2017/745 for medical devices; UMDNS code: 11-882.		
Materials:	Carboxylated nitrile butadiene rubber		
Product description:	<p>The nitrile glove is powder-free, latex-free and impervious to liquids.</p> <p>It offers excellent dexterity with high comfort and maximum dexterity (level 5/5). High elasticity with low effort. No hand fatigue.</p> <p>Ambidextrous wear, rolled edge, non-sterile, excellent tactile sensation.</p> <p>Textured surface all around the glove. Excellent grip, even when handling wet or moist objects.</p> <p>Excellent performance against chemicals, especially cytostatics and solvents, acids and bases, oxidizing agents, viruses, bacteria and fungi. Suitable for contact with foodstuff.</p> <p>EXTRA longer glove for even more safety against chemicals and EXTRA wearing comfort.</p>		
Pictograms and standards:	<div> <div>CAT 3</div> <div>ASTM 1671</div> <div>EN 21420</div> <div>EN 455-2</div> <div>AQL1,5</div> <div>ASTM D7160</div> <div>EN 455-1</div> <div>EN 455-3</div> <div> EN ISO 374-5:2016  VIRUS EN ISO 374-1:2016 TYPE C  </div> </div>		

Product performance:	<i>Testing against chemicals</i>		
	EN ISO 374-1:2016+A1:2018 and EN 16523-1:2015+A1:2018		
	Chemical	Level	EN ISO 374-4:2019 Degradation %

40% Sodium Hydroxide (K)	6	-3.9
30% Hydrogen Peroxide (P)	1	8.1
40% Acrylamide	6	NT
1% Ethidium Bromide	6	NT
NT = not tested		
Testing against infective agents accordingly EN ISO 374-5:2016 and ISO 16604:2004		
Protection against bacteria and fungi	Pass	
Protection against viruses	Pass	
Determination of resistance to penetration accordingly EN 374-2:2014		
Air leak test	Pass	
Water leak test	Pass	
Testing against the permeation of cytostatic drugs		
ASTM D6978-05(2019)		
Chemical	Minimum breakthrough detection time (Specimen 1/2/3) (Minutes)	Average steady state permeation rate (Specimen 1/2/3) (µg/cm²/Minute)
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	14.9 (16.9, 14.9, 14.3)	0.5 (0.6, 0.5, 0.5)
Cisplatin, 1 mg/ml (1,000 ppm)	>240	N/A
Cylophosphamid (Cytosan) 20.0 mg/ml (20,000 ppm)	>240	N/A
Dacarbazine 10.0 mg/ml (10,000 ppm)	>240	N/A
Doxorubicin HCl 2 mg/ml (2,000 ppm)	>240	N/A
Etoposide 20.0 mg/ml (20,000 ppm)	>240	N/A

	Fluorouracil 50.0 mg/ml (50,000 ppm)	>240	N/A
	Paclitaxel 6.0 mg/ml (6,000 ppm)	>240	N/A
	ThioTepa 10.0 mg/ml (10,000 ppm)	17.1 (17.2, 17.1, 18.5)	0.9 (0.9, 0.9, 0.8)
	<i>Dexterity and innocuousness of material testing</i>		
	EN ISO 21420:2020		
	Dexterity		Level 5
	pH Level		Pass
	Polycyclic Aromatic Hydrocarbons Content		Not detected - Pass
	<i>Food safety testing</i>		
	Resolution ResAP (2004) 4 and Regulation (EG)1935/2004		
	Overall migration		Pass
	Specific migration of primary aromatic amine		Pass
	Specific migration of nitrosamine and nitrosatable substances		Pass
	<i>Medical examination gloves testing</i>		
	Standard	Test	Result
	EN 455-1:2000	Freedom from holes	Pass
	EN 455-2:2015	Physical properties:	Pass
	Partial tests	Length: ≥ 240 mm	304 mm
		Width S: 80±10 mm	84 mm
		Width M: 95±10 mm	94 mm
		Width L: 110±10 mm	104 mm
		Width XL: ≥ 110 mm	113 mm
		Force at break: ≥ 6.0 N	9,25 N

		Force at break after challenge testing: ≥ 6.0 N	9,1 N
	EN 455-3:2015	Labelling	Pass
	Partial tests	Gloves shall not be dressed with talcum powder (magnesium silicate)	Pass
		Powder residues < 2mg	0,23 mg
	EN 455-4	Shelf life	Pass
Thickness:	Cuff: 0.09 mm; Palm: 0.1 mm; Finger: 0.14 mm (± 0.02 mm)		
Weight:	S: 5,5g, M: 6,0g, L: 6,8g, XL: 7,6g (± 0,4g)		
Other features:	<div>  <p>No intolerance reactions: This product does not contain natural latex.</p> </div> <div>  <p>Virus-tight according to ASTM 1671: Determination of the resistance of materials to the penetration of viral pathogens.</p> </div> <div>  <p>Non sterile.</p> </div>		
Color:	White		

Packaging and variants:	Packaging	Quantity
	Dispenser box	100 pcs.
	Shipping carton	1000 pcs.
GTIN:	879/W/S	4049825006688
	879/W/M	4049825006695
	879/W/L	4049825006701
	879/W/XL	4049825006718

Limitation:	<p>The gloves are designed for single use.</p> <p>The duration of protection in working use may differ from the breakthrough time determined according to EN 374, as it depends on the working conditions.</p>
Storage:	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 3 years (see labeling on the product).</p>
Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>