

Product Data Sheet



Product:	OP-Shorts plus	
Trademark:	SAFE®	
REF:	713	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.		
Materials:	▪ Nonwoven fabric: pure polypropylene ▪ Elastic: latex-free synthetic elastomer ▪ Sewing thread: 100 % polyester 17 - 18 tex		
Product description:	The OP-Shorts are made of soft non-woven fabric that is kind to the skin. It has half-length trouser legs, is equipped with elastic at the waist and has an examination slit with Velcro cover. The cut of the surgical shorts is generous and designed for double-sided use. Complies with the HACCP concept. Use: colonoscopy, gynecology, hematology.		
Dermatological compatibility:	No intolerance reactions.		
Dimensions:		Uni	XXL
	Length:	60 cm	67 cm
	Wide:	70 cm	80 cm
	Slot:	22,5 cm	25,5 cm
Colour:	Dark blue		
Packaging and variants:	Packaging	Quantity	
	Bag	10 pcs.	
	Shipping carton	100 pcs.	
GTIN:	713B/U	4049825001119	

	713B/XXL	4049825002130
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	