



**MADE IN
GERMANY**

**Top quality. Made in Germany
does not have to be expensive!**

DACH relies on a high degree of automation and intelligent product design for "Made in Germany".

DACH successful since 1996



**MADE IN
GERMANY**



Your strong partner since 1996 - more than a promise!

As a competent manufacturer of personal and medical protective equipment, our goal is to make the working world safer and to protect people as best as possible from damage to their health.

Our modern high-performance machines and our expertise offer the decisive advantage in quality and pricing. This is how we secure Germany as a production location in the long term.

You buy safe with us!

All DACH products are manufactured in strict compliance with European standards, inspected, packaged and documented. Our products offer high protection and excellent environmental compatibility. They are characterized by best quality, highest functionality and an excellent price-performance ratio.

In addition to our quality management system according to ISO 13485 and ISO 9001 and the requirements of legal regulations, our products are subjected to numerous additional tests. In the process, they have demonstrated the highest level of product safety and wearer comfort.





Modern, in-house test laboratory

Our in-house testing laboratory confirms our high quality standards. Regular monitoring of the goods and elaborate, continuous quality controls in the production process guarantee compliance with the highest quality standards.

So you can be sure that our products more than fulfill their protective functions.

Worldwide competitive due to high degree of automatization

Company since 1996 with high innovative power & professional expertise



Internationally competitive prices thanks to high degree of automatization



Manufacturing in Rastatt, Germany: Strengthening the production location



Regionality increases sustainability and reduces the CO2 footprint



WE KEEP YOU **SAFE**[®]
Since 1996

neolution[®] Air

Highest filter performance, feather-light and extremely air-permeable - for a cool feeling when worn.



**MADE IN
GERMANY**



The high-end FFP mask is a perfect combination of protective effect, wearing comfort and beautiful design. It meets the increased demands today placed on protective masks.

Easy breathing, maximum protection through extremely high filtration – the neolution® Air performs

**Passed Dolomit-Test
Certified
biocompatibility
according to EN 10993**

Integrated nose clip: Perfect fit, no direct contact with metal

Intelligent 3-zone design ensures a tight fit on the face

Very easy breathing thanks to special fabric

Ultralight: Only 3.5 grams with the highest level of protection



Suitable for a long time of usage thanks to generous breathing volume


Chic flat design with colorful accents for style-conscious wearers
Biological risk assessment and compatibility according to EN 10993




REF	Protection Class	Packed	Carton
255	FFP2 NR D, EN 10993	50 Pcs./Box	300 Pcs.
256	FFP3 NR D, EN 10993	50 Pcs./Box	300 Pcs.
255-3	FFP2 NR D, EN 10993	3 Pcs./Bag	270 Pcs.

REF	Protection Class	Packed	Carton
253	FFP2 NR D, EN 10993	50 Pcs./Box	300 Pcs.
257	FFP3 NR D, EN 10993	50 Pcs./Box	300 Pcs..
253-3	FFP2 NR D, EN 10993	50 Pcs./Box	270 Pcs.

Average inhalation resistance of DACH respirators - according to DIN EN 149:2001+A1:2009



Average NaCl permeability of DACH respirators - according to DIN EN 149:2001+A1:2009



ubiQ[®]

Combines high quality, comfort and great economy for ubiquitous, professional use

Passed Dolomit-Test
FFP2 NR D & FFP3 NR D


 **MADE IN GERMANY**

The large filter surfacefl che ensures excellent particle distribution and thus provides maximum protection and high wearing comfort even over longer durations of wear.

- Tested according to DIN EN 149 and DIN EN 14683.
- Ergonomic fit and soft material guarantee a tight fit for any face shape.
- High filter performance with low breathing resistance.
- 2-piece design for easy handling.

REF	Protection Class	Packed	Carton
272	FFP2 NR D, Type IIR, EN 10993	50 Pcs./Box	900 Pcs.
273	FFP3 NR D, Type IIR, EN 10993	50 Pcs./Box	900 Pcs.

Average NaCl permeability of DACH respirators - according to DIN EN 149:2001+A1:2009



Respirator and medical face mask in one. FFP2 or FFP3.
Whenever self-protection and external protection have to be combined, whenever the optimum protective mask must be selected: The ubiQ[®] is definitely the best choice.

Model D Loop

Medical face mask with ear loops

**MADE IN
GERMANY**



Made with high performance machines in Germany.

- Liquid repellent
- Low breathing resistance
- Skin friendly hypoallergenic material
- Best price-performance ratio

REF	Protection Class	Packed	Carton
216	Typ IIR, EN 10993	50 Pcs./Box	1.100 Pcs.
217	Typ IIR, EN 10993	5 Pcs./Bag	300 Pcs.

Model D Loop in vielen, aktuellen Farben erhältlich.



REF: 216: black, purple, blue, rosa

REF: 217: blue, orange, red, rosa, beige

Customized packaging upon request

Low price and good quality - this is what characterizes our Model D Loop. The high degree of automation of our production ensures our international competitiveness.

WE KEEP YOU **SAFE**[®]
SINCE 1996

DACH
Schutzbekleidung



MADE IN
GERMANY



RESPIRATORY



PROTECTIVE CLOTHING



HAND PROTECTION



HEADGEAR



DACH Schutzbekleidung GmbH & Co. KG

Rotackerstraße 21
76437 Rastatt, Germany

Phone: +49 7222 15989-0
Fax: +49 7222 15989-22

Mail: info@dach-germany.de
Web: www.dach-germany.de