

Product:	ClipCap	
Trademark:	ecoComfort®	
REF:	102	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.	
Materials:	<ul style="list-style-type: none"> ▪ Nonwoven fabric: pure polypropylene ▪ Elastic: latex-free synthetic elastomer 	
Product Description:	<p>The cap is made of soft nonwoven that has been folded into a strip. Pulling the folds apart creates a round cap.</p> <p>The outer finish is provided with a welded-in elastic.</p> <p>Conforms to the HACCP concept.</p>	
Dermatological compatibility:	No intolerance reactions.	
Colours:	Blue, White	
Packaging and variants:	Packaging	Quantity
	Dispenser box	100 pcs.
	Shipping carton	1.000 pcs.
GTIN:	102B	4049825001638
	102W	4049825001584
Storage:	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 5 years. (see labeling on the product)</p>	
Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	