


Product:	Bouffant Cap	
Trademark:	SAFE®	
REF:	113	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I Personal Protective Equipment accordingly Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.	
Materials:	<ul style="list-style-type: none"> ▪ Hood: pure polypropylene ▪ Elastic band: latex-free, synthetic elastomer ▪ Sewing thread: 100 % Polyester 	
Product Description:	<p>The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam.</p> <p>In particular, protection is provided against tangling of hair in the immediate working environment of the wearer.</p> <p>Conforms to the HACCP concept.</p>	
Dermatological tolerance:	No intolerance reactions.	
Standard size:	53 cm	
Colours:	Blue, White, Green	
Packaging and Variants:	Quantity	Packaging
	200 pcs.	Bag
	1.000 pcs.	Shipping carton
GTIN:	113B	4049825000303
	113W	4049825000327
	113G	4049825000310
Storage:	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).</p>	



Environmental sustainability and disposal:	<p>In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>
---------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------