

Product:	Bouffant Cap	
Trademark:	SAFE®	
REF:	113	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	▪ CAT I Personal Protective Equipment accordingly Regulation (EU) 2016/425. ▪ Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223. ▪ Helmke Drum Test accordingly IEST-PR-CC003.4: Category II We recommend the use in clean rooms accordingly EN ISO 14644-1: ISO class 6 to 9.	
Materials:	▪ Hood: pure polypropylene ▪ Elastic band: latex-free, synthetic elastomer ▪ Sewing thread: 100 % Polyester	
Product Description:	The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam. In particular, protection is provided against tangling of hair in the immediate working environment of the wearer. Conforms to the HACCP concept.	
Dermatological tolerance:	No intolerance reactions.	
Standard size:	53 cm	
Colours:	Blue, White, Green	
Packaging and Variants:	Quantity	Packaging
	200 pcs.	Bag
	1.000 pcs.	Shipping carton
GTIN:	113B	4049825000303
	113W	4049825000327

	113G	4049825000310
Storage:	Store in a dry place without direct sunlight in the original packaging (see packaging). If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).	
Environmental sustainability and disposal:	In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	