


<b>Product:</b>	PP-Shoe Covers	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	931	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Nonwoven fabric: pure polypropylene</li> <li>▪ Elastic: high quality polyurethane, latex free</li> <li>▪ Sewing thread: 100 % polyester</li> </ul>
<b>Product description:</b>	<p>The overshoes are made of strong non-woven fabric.</p> <p>The upper opening of the shoes is formed by an elastic band bordered with an overlock seam.</p> <p>The sole of the foot is rubberized.</p> <p>Conforms to the HACCP concept.</p>
<b>Dermatological compatibility:</b>	No intolerance reactions
<b>Size:</b>	Length: 44 cm, width: 16 cm.
<b>Color:</b>	Blue, White

Packaging and variants:	Packaging	Quantity
	Bag, 10 pcs rolled up	100 pcs.
	Shipping carton	1.000 pcs.
<b>GTIN:</b>	931/B	4049825011743
	931/W	4049825011736
<b>Storage:</b>	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>	
<b>Environmental sustainability and disposal:</b>	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	