

Product:	PP-Shoe Covers	
Trademark:	SAFE®	
REF:	931	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.	
Materials:	▪ Nonwoven fabric: pure polypropylene ▪ Elastic: high quality polyurethane, latex free ▪ Sewing thread: 100 % polyester	
Product description:	The overshoes are made of strong non-woven fabric. The upper opening of the shoes is formed by an elastic band bordered with an overlock seam. The sole of the foot is rubberized. Conforms to the HACCP concept.	
Dermatological compatibility:	No intolerance reactions	
Size:	Length: 44 cm, width: 16 cm.	
Color:	Blue, White	
Packaging and variants:	Packaging	Quantity
	Bag, 10 pcs rolled up	100 pcs.
	Shipping carton	1.000 pcs.
GTIN:	931/B	4049825011743
	931/W	4049825011736
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	