

<b>Product:</b>	Bouffant cap	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	115	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I Personal Protective Equipment accordingly Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.	
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Hood: pure polypropylene</li> <li>▪ Elastic band: latex-free, synthetic elastomer</li> <li>▪ Sewing thread: 100 % Polyester</li> </ul>	
<b>Product description:</b>	<p>The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam.</p> <p>In particular, protection is provided against tangling of hair in the immediate working environment of the wearer.</p> <p>Conforms to the HACCP concept.</p>	
<b>Dermatological tolerance:</b>	No intolerance reactions.	
<b>Size:</b>	61 cm	
<b>Colours:</b>	Blue, White, Green	
<b>Packaging and Variants:</b>	<b>Packaging</b>	<b>Quantity</b>
	Box	200 pcs.
	Shipping carton	1.000 pcs.
<b>GTIN:</b>	115B	4049825000365
	115W	4049825000389
	115G	4049825000372
<b>Storage:</b>	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).</p>	



<b>Environmental sustainability and disposal:</b>	<p>In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>
---	---