

Product:	Bouffant cap	
Trademark:	SAFE®	
REF:	115	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I Personal Protective Equipment accordingly Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.	
Materials:	▪ Hood: pure polypropylene ▪ Elastic band: latex-free, synthetic elastomer ▪ Sewing thread: 100 % Polyester	
Product description:	The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam. In particular, protection is provided against tangling of hair in the immediate working environment of the wearer. Conforms to the HACCP concept.	
Dermatological tolerance:	No intolerance reactions.	
Size:	61 cm	
Colours:	Blue, White, Green	
Packaging and Variants:	Packaging	Quantity
	Box	200 pcs.
	Shipping carton	1.000 pcs.
GTIN:	115B	4049825000365
	115W	4049825000389
	115G	4049825000372
Storage:	Store in a dry place without direct sunlight in the original packaging (see packaging). If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).	



Environmental sustainability and disposal:	In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.
---	---