


Product:	Disposable sheet	
Trademark:	ecoComfort®	
REF:	3008	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Materials:	<ul style="list-style-type: none"> ▪ Pulp ▪ Polyethylene ▪ Nylon thread 	
Product description:	The carrying sheet consists of a combination of absorbent cellulose and a liquid-impermeable film. The 48-fold thread reinforcement guarantees extra high tear strength.	
Size:	80 x 210 cm	
Packaging and variants:	Packaging	Quantity
	Bag	25 pcs.
	Shipping carton	100 pcs.
GTIN:	3008	4049825002970
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	