


Product:	MediCap AstroEco	
Trademark:	SAFE®	
REF:	122G	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	Class I of Directive EU 2017/745 on medical devices ; UMDNS code: 15-223.	
Materials:	<ul style="list-style-type: none"> ▪ Hood: pure polypropylene nonwoven, antistatic treated ▪ Elastic band: latex-free, synthetic elastomer ▪ Sewing thread: 100 % Polyester 	
Product Description:	<p>The hood in astronaut shape consists of a head top and a round part.</p> <p>The two parts are joined together with an overlock seam.</p> <p>The all-around part can also be used for tying.</p> <p>The front side of the hood is provided with a multi-layered, absorbent sweatband.</p>	
Dermatological compatibility:	No intolerance reactions.	
Colour:	Green	
Packaging and Variants:	Packaging	Quantity
	Bag	100 pcs.
	Shipping carton	1.000 pcs.
GTIN:	122G	4049825000464
Storage:	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).</p>	
Environmental sustainability and disposal:	<p>In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	