


<b>Product:</b>	MediCap AstroEco, with sweat band	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	122G	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	Class I of Directive EU 2017/745 on medical devices ; UMDNS code: 15-223.	
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Hood: pure polypropylene nonwoven, antistatic treated</li> <li>▪ Elastic band: latex-free, synthetic elastomer</li> <li>▪ Sewing thread: 100 % Polyester</li> </ul>	
<b>Product Description:</b>	<p>The hood in astronaut shape consists of a head top and a round part.</p> <p>The two parts are joined together with an overlock seam.</p> <p>The all-around part can also be used for tying.</p> <p>The front side of the hood is provided with a multi-layered, absorbent sweatband.</p>	
<b>Dermatological compatibility:</b>	No intolerance reactions.	
<b>Colour:</b>	Green	
<b>Packaging and Variants:</b>	<b>Packaging</b>	<b>Quantity</b>
	Bag	100 pcs.
	Shipping carton	1.000 pcs.
<b>GTIN:</b>	122G	4049825000464
<b>Storage:</b>	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).</p>	
<b>Environmental sustainability and disposal:</b>	<p>In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	