

Product:	MediComfort Excel	
Trademark:	SAFE®	
REF:	221G	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	Class I of Regulation (EU) 2017/745 on medical devices; UMDNS Code: 12-458			
Protection class:	Type II R (liquid resistant), tested according to EN 14683:2019 for surgical masks.			
Materials:	▪ Outside: liquid-repellent polypropylene fleece ▪ Inner side: especially soft and skin-friendly, allergen-free polypropylene fleece ▪ Filter medium: microfiber made of polypropylene ▪ Straps: skin-friendly, allergy-free polypropylene fleece ▪ Noseband: PE coated wire			
Product Description:	Latex-free, high-quality, anatomically shaped face mask with easy-to-mold nosepiece. The inside and outside are made of physiologically harmless unplasticized polypropylene fleece. The filter medium is a special microfiber fleece manufactured in the meltblown process, welded together in the ultrasonic process.			
Product performance:	Product performance according to EN 14683:2019+AC:2019			
	Testing	Performance class	Requirement	Test result
	Bacterial filtration efficiency (BFE), (%))	Type IIR	≥ 98	99,9
	Pressure difference (Pa/cm²)		< 60	56,5
	Spray resistance pressure (kPa)		≥ 16,0	16
	Microbial cleanliness (CFU/g)		≤ 30	≤ 25,6
Dermatological compatibility:	No intolerance reactions			
Colour:	Green striped			
Packaging and variants:	Packaging		Quantity	

	Dispenser box	50 pcs.
	Shipping carton	1.000 pcs.
GTIN:	221G	4049825000198
Use:	All materials used are free of irritants, odorless and suitable for medical purposes. The mask is suitable for permanent use in clean rooms, e.g. operating rooms, pharmaceuticals, food and electronics, and provides excellent protection.	
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 5 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	
Indication:	Non sterile. 	