


# Product Data Sheet



<b>Product:</b>	Bouffant Cap	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	111	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.	
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Cap: 100% polypropylene non-woven fabric</li> <li>▪ Elastic: polyurethane elastomere (free from natural latex)</li> <li>▪ Sewing thread: 100 % polyester</li> </ul>	
<b>Product description:</b>	<p>The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam.</p> <p>In particular, protection is provided against tangling of hair in the immediate working environment of the wearer.</p> <p>Conforms to the HACCP concept.</p>	
<b>Dermatological compatibility:</b>	No intolerance reactions.	
<b>Size:</b>	53 cm	
<b>Colours:</b>	Blue, White, Green	
<b>Packaging and variants:</b>	<b>Packaging</b>	<b>Quantity</b>
	Bag	100 pcs.
	Shipping carton	1.000 pcs.
<b>GTIN:</b>	111B	4049825000334
	111W	4049825000358
	111G	4049825000341
<b>Storage:</b>	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>	

**Environmental compatibility and disposal:**

In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.

A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.