

|                      |                                     |   |
|----------------------|-------------------------------------|---|
| <b>Product:</b>      | Bouffant Cap                        |  |
| <b>Trademark:</b>    | SAFE®                               |   |
| <b>REF:</b>          | 111                                 |   |
| <b>Manufacturer:</b> | DACH Schutzbekleidung GmbH & Co. KG |   |

|                                      |   |                 |
|--------------------------------------|---|-----------------|
| <b>Classification:</b>               | <ul style="list-style-type: none"> <li>CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.</li> <li>Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 15-223.</li> <li>Helmke Drum Test accordingly IEST-PR-CC003.4: Category II</li> </ul> <p>We recommend the use in clean rooms accordingly EN ISO 14644-1: ISO class 6 to 9.</p> |                 |
| <b>Materials:</b>                    | <ul style="list-style-type: none"> <li>Cap: 100% polypropylene non-woven fabric</li> <li>Elastic: polyurethane elastomere (free from natural latex)</li> <li>Sewing thread: 100 % polyester</li> </ul>  |                 |
| <b>Product description:</b>          | <p>The cap is made of soft nonwoven in the form of a bouffant cap. The outer finish of the cap is an elastic seamed by an overlock seam.</p> <p>In particular, protection is provided against tangling of hair in the immediate working environment of the wearer.</p> <p>Conforms to the HACCP concept.</p>  |                 |
| <b>Dermatological compatibility:</b> | No intolerance reactions.   |                 |
| <b>Size:</b>                         | 53 cm   |                 |
| <b>Colours:</b>                      | Blue, White, Green  |                 |
| <b>Packaging and variants:</b>       | <b>Packaging</b>  | <b>Quantity</b> |
|                                      | Bag   | 100 pcs.        |
|                                      | Shipping carton   | 1.000 pcs.      |
| <b>GTIN:</b>                         | 111B  | 4049825000334   |
|                                      | 111W  | 4049825000358   |
|                                      | 111G  | 4049825000341   |

|  |  |
|--|--|
| <b>Storage:</b>                                  | <p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>   |
| <b>Environmental compatibility and disposal:</b> | <p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p> |