

Product:	Visitor's Set	
REF:	9015	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Product Description:	The visitor's set includes a disposable gown and a disposable bouffant cap, hygienically sealed in a bag. Complies with the HACCP concept. Available in white or red. The protective set consists of the following components:	
Components:	1 Gown:	The inspection gown is made of soft nonwoven fabric, with an opening in the front, closed with a Velcro fastener. The sleeve ending is finished with an elastic (free from natural rubber) enclosed by an overlock seam. PPE Cat. I
	1 Hood:	The bouffant cap is made of soft nonwoven fabric. The outer edge at the back of the hood is finished with an elastic (free from natural latex) enclosed by an overlock seam. The cap has a sewn-in peak at the front. PPE Cat. I

Packaging and Variants:	REF	Color	Packaging	Quantity	GTIN
	9015/W	White	Bag	1 pc.	4049825000235
	9015/W	White	Shipping Carton	25 pcs.	4049825500223
	9015/R	Red	Bag	1 pc.	4049825004400
	9015/R	Red	Shipping Carton	25 pcs.	4049825504399
Euro Pallet: 1600 pcs.					
Storage:	Store dry in original packaging, away from direct sunlight, see also packaging. Maximum shelf life: 10 years				



Disposal:	The product can be thermally recycled or disposed of in landfills without releasing toxic substances. For contaminated products, the type and extent of contamination determine disposal, and applicable laws and regulations must be followed.
------------------	--