



<b>Product:</b>	Patient Transport+ protection kit	
<b>REF:</b>	9017	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Product description:</b>	<p>The Patient Transport+ Protection Set was specially developed for use in clinical transport. It offers the user extended protection against viruses, bacteria and spray of physical and chemical fluids.</p> <p>The protective set consists of the following components:</p>	
<b>Components:</b>	1 Protective suit:	 <p>1 Full body protective suit KAT III, type 4b+5b+6b, antistatic according to EN 1149-1</p>
	1 Respirator mask:	"High-risk mask", protection class FFP3 NR D, without exhalation valve, for protection against airborne infectious agents
	1 pair of overshoes:	PE overshoes, liquid-tight, made of stable PE film
<b>Storage:</b>	<p>Store without direct sunlight, dry in original packaging, see also packaging.</p> <p>Maximum shelf life: 3 years</p>	
<b>Disposal:</b>	<p>The product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p> <p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations must be observed.</p>	
<b>Warnings:</b>	Before using the set, read the inside instructions for use of the respective components.	