

Product:	OP-Shorts	
Trademark:	SAFE®	
REF:	712	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.		
Materials:	<ul style="list-style-type: none"> ▪ Nonwoven fabric: 100% polypropylene ▪ Elastic: polyurethane - elastomers (free from natural latex) ▪ Sewing thread: 100 % polyester 		
Product description:	<p>The OP-Shorts are made of soft non-woven fabric that is kind to the skin.</p> <p>It has half-length trouser legs, is equipped with elastic at the waist and has an examination slit with overlapping fabric.</p> <p>The cut of the surgical shorts is generous and designed for double-sided use. Conforms to the HACCP concept.</p> <p>Use: colonoscopy, gynecology, hematology.</p>		
Dermatological compatibility:	No intolerance reactions.		
Dimensions:		Uni	XXL
	Length:	65 +-2 cm	65 +-2 cm
	Wide:	60 +-2 cm	75 +-2 cm
	Slot:	37 +-0,5 cm	37 +-0,5 cm
Color:	Dark blue		



Packaging and variants:	Packaging	Quantity
	Bag	10 pcs.
	Shipping carton	100 pcs.
GTIN:	712B/U	4049825001119
	712B/XXL	4049825002130
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	