


Product:	PP-Shoe Covers	
Trademark:	SAFE®	
REF:	925W	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I of Regulation (EU) 2016/425 for personal protective equipment (PPE).
Materials:	<ul style="list-style-type: none"> Shoe: pure polypropylene, antistatic treated Sole: polyethylene cover Elastic band: Natural rubber
Product description:	The shoe cover is made of resilient non-woven fabric. The lower half is coated with a liquid-repellent PE coating. All seams have been firmly welded using ultrasonic technology.
Size:	Length: 41 cm, width: 15 cm.
Color:	Upper white / sole blue.

Packaging and variants:	Packaging	Quantity
	Bag, 5 pcs rolled up	100 pcs.
	Shipping carton	500 pcs.
GTIN:	925W	4049825002307
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	