


<b>Product:</b>	PP-Shoe Covers	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	926W	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425.
<b>Material:</b>	<ul style="list-style-type: none"> <li>Shoe: pure polypropylene, antistatic treated</li> <li>Sole: polyethylene cover</li> <li>Elastic band: Natural rubber</li> </ul>
<b>Product description:</b>	The shoe cover is made of resilient non-woven fabric. The lower half is coated with a liquid-repellent PE coating. All seams have been firmly welded using ultrasonic technology.
<b>Size:</b>	Length: 41 cm, width: 15 cm.
<b>Color:</b>	Upper white / sole white

Packaging and variants:	Packaging	Quantity
	Bag, 5 pcs rolled up	50 pcs.
	Shipping carton	500 pcs.
<b>GTIN:</b>	926W	4049825003397
<b>Storage:</b>	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
<b>Environmental sustainability and disposal:</b>	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	