

Product:	Biohazard disposal bag	
REF:	017B	

Product description:	Disposal bag according to the quality specifications of TRBA 250 (Technical Rules for Biological Agents 250) and DGUV Information 203-084 (Information 203-084 of the German statutory accident insurance). The bag is not autoclavable.	
Material:	LDPE	
Dimensions:	Thickness: 0,085 mm Width: 550 mm Length: 1000 mm Volume: ca. 70 l	
Colour:	Yellow with black print.	
Packaging and Variants:	Quantity	Packaging
	20 pcs.	Roll
	200 pcs.	Shipping Carton
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging.	