

|                      |                                     |                                                                                     |
|----------------------|-------------------------------------|-------------------------------------------------------------------------------------|
| <b>Product:</b>      | PE-Sleeves                          |  |
| <b>Trademark:</b>    | SAFE®                               |                                                                                     |
| <b>REF:</b>          | 651                                 |                                                                                     |
| <b>Manufacturer:</b> | DACH Schutzbekleidung GmbH & Co. KG |                                                                                     |

|                                      |                                                                                                                                                                                                                                                                                             |                 |
|--------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| <b>Classification:</b>               | <ul style="list-style-type: none"> <li>▪ CAT I of Regulation (EU) 2016/425 for personal protective equipment (PPE).</li> <li>▪ Helmke Drum Test accordingly IEST-PR-CC003.4: Category I</li> </ul> <p>We recommend the use in clean rooms accordingly EN ISO 14644-1: ISO class 5 to 9.</p> |                 |
| <b>Materials:</b>                    | <ul style="list-style-type: none"> <li>▪ Sleeves: 100 % polyethylene</li> <li>▪ Elastic band: Natural rubber</li> </ul>                                                                                                                                                                     |                 |
| <b>Product description:</b>          | <p>The sleeves are made of a liquid impermeable and recyclable film, tubular with elastic band welded together at both ends.</p> <p>Complies with the HACCP concept.</p>                                                                                                                    |                 |
| <b>Dermatological compatibility:</b> | No intolerance reactions                                                                                                                                                                                                                                                                    |                 |
| <b>Size:</b>                         | Length: 40 cm, Wide: 20 cm.                                                                                                                                                                                                                                                                 |                 |
| <b>Colour:</b>                       | Blue, White                                                                                                                                                                                                                                                                                 |                 |
| <b>Packaging and variants:</b>       | <b>Packaging</b>                                                                                                                                                                                                                                                                            | <b>Quantity</b> |
|                                      | Bag, 20 pcs rolled up                                                                                                                                                                                                                                                                       | 100 pcs.        |
|                                      | Shipping carton                                                                                                                                                                                                                                                                             | 2000 pcs.       |
| <b>GTIN:</b>                         | 651B                                                                                                                                                                                                                                                                                        | 4049825001096   |
|                                      | 651W                                                                                                                                                                                                                                                                                        | 4049825001102   |
| <b>Storage:</b>                      | <p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>                                                                                    |                 |

**Environmental  
sustainability and  
disposal:**

In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.

A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.