


Product:	Medical Face Mask	
REF:	204	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 12-458
Protection class:	Type IIR, tested according to EN 14683:2019 for medical face masks.
Materials:	<ul style="list-style-type: none"> ▪ Outside: polypropylene ▪ Inner side: skin-friendly, allergy-free polypropylene ▪ Filter medium: polypropylene ▪ Straps: soft polyester and elastane ▪ Noseband: plastic coated wire
Product description:	<p>Latex-free, allergy-free, anatomically shaped face mask with easy-to-mold nosepiece.</p> <p>The inside and outside are physiologically harmless and plasticizer-free.</p> <p>Ultrasonically welded and with elastic ear loops.</p>
Dermatological compatibility:	No intolerance reactions.

Product performance:	Product performance according to EN 14683:2019+AC:2019			
	Testing	Performance class	Requirement	Test result
	Bacterial filtration efficiency (BFE), (%)	Typ IIR	≥ 98	99,9
	Pressure difference (Pa/cm²)		< 60	36,6
	Spray resistance pressure (kPa)		≥ 16,0	16
Colour:	Black			
Packaging and variants:	Packaging		Quantity	
	Dispenser box		50 pcs.	

	Shipping carton	2.000 pcs.
GTIN:	204	4049825006459
Use:	All materials used are free of irritants, odorless and suitable for medical purposes. By wearing the mask, the wearer protects his environment.	
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 5 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	
Indication:	Non sterile 