

Product Data Sheet



Product:	Disposable warm-up jacket	
REF:	753	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	▪ Class I, based on Medical Device Regulation (EU) 2017/745 MDR, Annex VIII ▪ Helmke Drum Test accordingly IEST-PR-CC003.4: Category II We recommend the use in clean rooms accordingly EN ISO 14644-1: ISO class 7 to 9.		
Materials:	▪ Nonwoven fabric: pure polypropylene ▪ Knitted cuff: 100 % polyester ▪ Collar: 100 % cotton ▪ Sewing thread: 100 % polyester ▪ Press studs: polyoxymethylene		
Product description:	The disposable warm-up jacket made of robust, low linting SMS trilaminate excels wherever personal clothing cannot be used due to increased hygiene requirements. The warming material is breathable and water-repellent. The jacket does not "stick" to the body due to a special antistatic treatment. Long sleeves with knitted cuffs and the soft knitted collar at the neck provide high wearing comfort. The jacket has three pockets. The product is not sterile.		
Dimensions:	Size S: Length: 64 cm, Width: 48 cm Size M: Length: 72 cm, Width: 55 cm Size L: Length: 80 cm, Width: 62 cm Size XL: Length: 88 cm, Width: 69 cm Size XXL: Length: 96 cm, Width: 76 cm Size XXXL: Length: 104 cm, Width: 83 cm The dimensions are subject to a tolerance of +/- 2 cm.		
Colour:	Blue, Green		
Packaging and variants:	Packaging	Quantity	
	Bag	10 pcs.	
	Shipping carton	50 pcs.	
GTIN:	REF	Primary packaging (Bag)	Secondary packaging (Shipping carton)
	753/B/S	4049825008828	4049825508816
	753/B/M	4049825008835	4049825508823

	753/B/L	4049825008842	4049825508830
	753/B/XL	4049825008859	4049825508847
	753/B/XXL	4049825008866	4049825508854
	753/B/XXXL	4049825008873	4049825508861
	753/G/S	4049825008880	4049825508878
	753/G/M	4049825008897	4049825508885
	753/G/L	4049825008903	4049825508892
	753/G/XL	4049825008910	4049825508908
	753/G/XXL	4049825008927	4049825508915
	753/G/XXXL	4049825008934	4049825508922
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 5 years. (see labeling on the product)		
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.		