



<b>Product:</b>	Protection kit Pro+	
<b>REF:</b>	9003	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Product description:</b>	<p>This set of personal protective equipment provides the greatest protection in both professional and industrial environments.</p> <p>The protective set consists of the following components:</p>	
<b>Components:</b>	1 Protective suit:	<p>KAT3 EN1149-1 TYP 5 TYP 6</p>  <p>1 Overall, PSA Kat. 3, Typ 5+6</p>
	1 Respirator mask:	Protection class FFP3 D, with exhalation valve, additional dolomite test, also for industrial applications
	2 pairs of examination gloves:	PPE cat. 3 according to EN 374 (in case of lethal risks), liquid- and virus-proof, great wearing comfort
	1 pair boot protectors	PPE cat. 1, liquid-tight, overlapped with coverall
	1 Safety goggles	PPE cat. 2, enclosing, indirect ventilation, EN 166
	1 Disposal bag	After use, all parts can be collected in the bag and taken for disposal. Disposal must be done professionally according to the type of contamination.

Packaging and variants:	Packaging	Quantity
	Bag	1 pc.
	Shipping carton	25 pcs.

<b>GTIN:</b>	9003/M	4049825002697
	9003/L	4049825001232
	9003/XL	4049825001249
	9003/XXL	4049825001256
<b>Storage:</b>	Store without direct sunlight, dry in original packaging, see also packaging. Maximum shelf life: 3 years	
<b>Disposal:</b>	The product can be thermally recycled or disposed of in landfills without releasing toxic substances. In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations must be observed.	
<b>Warnings:</b>	Before using the set, read the inside instructions for use of the respective components.	

## **Information: Protection against Ebola**

Due to the strong demand, DACH Schutzbeleidung GmbH & Co. KG has developed three different Ebola protection sets based on information from the Robert Koch Institute (RKI), which are tailored to the needs in different situations when dealing with sick people. These situations are (the word patient may also stand for suspected disease in the following):

- I. Keeping a minimum distance of 1.20 m from the patient, but entering an area that may be contaminated. Splashing of blood, secretions, or excretions is not expected.
- II. the minimum distance of 1.20 m to the patient is not maintained, but the patient is not treated directly. The protective measures serve to protect against contact and the splashing of blood, secretions or excrement is not to be expected.
- III. direct treatment of the patient is necessary or splashing of blood, secretions or excrements is to be expected.

**We recommend the Profi+ protective set for situation II. It is the customer's responsibility to check the professional applicability of the products under the given circumstances.**