

# Product Data Sheet



<b>Product:</b>	Surgical Gown Easy	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	443	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I of Regulation (EU) 2016/425 for personal protective equipment (PPE). Class I of Regulation (EU) 2017/745 for medical devices.
<b>Protection class:</b>	DIN EN 13795:2011+A1:2013 Standard Performance
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Sheath: SMS</li> <li>▪ Outside: Liquid repellent nonwoven fabric</li> <li>▪ Inner side: Skin-friendly nonwoven fabric</li> <li>▪ Filter layer: Microfiber fleece</li> <li>▪ Knitted cuff: 100 % polyester</li> </ul> <p>Free from PVC (polyvinyl chloride), DEHP(2-ethylhexyl phthalate), latex</p> <p><b>This product is not sterile.</b></p>
<b>Size:</b>	Length: 43 cm, width: 20 cm.

<b>Packaging and variants:</b>	<b>Packaging</b>	<b>Quantity</b>
		Bag
	Shipping carton	50 pcs.
<b>GTIN:</b>	443/G/L	4049825003489

<b>Storage:</b>	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)
<b>Environmental sustainability and disposal:</b>	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.