

<b>Product:</b>	IsoCoat Zytex® 360	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	425Y	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 11-897.	
<b>Materials:</b>	<ul style="list-style-type: none"> <li>▪ Gown: Zytex®</li> <li>▪ Sleeve cuffs: knitted cuff made of 100 % polyester</li> <li>▪ Sewing thread: 100 % polyester</li> </ul>	
<b>Product description:</b>	<p>The isolation gown is made of Zytex® material, which has been tested against the penetration of blood and viruses according to EN 14126 and against the permeation of cytostatics according to EN ISO 6529.</p> <p>Zytex® was specially developed for handling cytostatic drugs and has a high impermeability to cytostatic drugs.</p> <p>Opening on the back, two tie belts and with knitted cuffs.</p>	
<b>Dermatological compatibility:</b>	No intolerance reactions.	
<b>Size/Length:</b>	130cm	
<b>Colour:</b>	Yellow	

<b>Packaging and variants:</b>	<b>Packaging</b>	<b>Quantity</b>
	Bag	10 pcs.
	Shipping carton	50 pcs.
<b>GTIN:</b>	425Y	4049825002444
<b>Use:</b>	Infection prevention and protection against MRSA or noroviruses	

<b>Storage:</b>	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>
<b>Environmental sustainability and disposal:</b>	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>