

Product:	IsoCoat Zytex® 360	
Trademark:	SAFE®	
REF:	425Y	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	<ul style="list-style-type: none"> ▪ CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. ▪ Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 11-897. ▪ Helmke Drum Test accordingly IEST-PR-CC003.4: Category II <p>We recommend the use in clean rooms accordingly EN ISO 14644-1: ISO class 6 to 9.</p>	
Materials:	<ul style="list-style-type: none"> ▪ Gown: Zytex® ▪ Sleeve cuffs: knitted cuff made of polyester ▪ Sewing thread: polyester 	
Product description:	<p>The isolation gown is made of Zytex® material, which has been tested against the penetration of blood and viruses according to EN 14126 and against the permeation of cytostatics according to EN ISO 6529.</p> <p>Zytex® was specially developed for handling cytostatic drugs and has a high impermeability to cytostatic drugs.</p> <p>Opening on the back, two tie belts and with knitted cuffs.</p>	
Dermatological compatibility:	No intolerance reactions.	
Size/Length:	130cm	
Colour:	Yellow	
Packaging and variants:	Packaging	Quantity
	Bag	10 pcs.
	Shipping carton	50 pcs.
GTIN:	425Y	4049825002444
Use:	Infection prevention and protection against MRSA or noroviruses	

Storage:	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>
Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>