

Product:	EcoCoat pro	
Trademark:	ecoComfort®	
REF:	411	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 11-897.	
Materials:	<ul style="list-style-type: none"> ▪ nonwoven fabric: polypropylene ▪ Sleeve cuffs: polyester ▪ Sewing thread: polyester 	
Product description:	<p>The Eco-Coat pro is made of soft non-woven fabric, with an opening on the back, two tie belts and knitted cuffs.</p> <p>Complies with the HACCP concept.</p>	
Dermatological compatibility:	No intolerance reactions.	
Length:	Size M:115 cm, size L: 136 cm, size XL: 150 cm	
Colour:	Green	
Packaging and variants:	Packaging	Quantity
	Bag	10 pcs.
	Shipping carton	100 pcs.
GTIN:	411/G/M	4049825002703
	411/G/L	4049825003915
Storage:	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>	
Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	