


Product:	StandSafe	
Trademark:	SAFE®	
REF:	976	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I of Regulation (EU) 2016/425 for personal protective equipment (PPE).
Materials:	<ul style="list-style-type: none"> ▪ SMS: - Outside: Liquid repellent polypropylene fleece - Inside: Skin-friendly, allergen-free polypropylene fleece - Filter media: Microfiber produced in meltblown process. ▪ Sole: ethylene vinyl acetate (copolymer), slip-resistant ▪ Elastic: Latex-free, synthetic elastomer
Product description:	<p>The StandSafe consists of a non-slip sole and a particle-filtering, liquid-repellent material made of SMS.</p> <p>The opening of the overshoe is bordered with an elastic band.</p> <p>The upper material is electrostatically dissipative according to EN 1149-1.</p>
Dermatological compatibility:	No intolerance reactions
Size:	<p>Length: 36 cm, height: 18cm</p> <p>Suitable for shoe size 37 to 46.</p>
Color:	White

Packaging and variants:	Packaging	Quantity
	Bag	10 Pair
	Shipping carton	200 Pair
GTIN:	976	4049825004622
Storage:	<p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p>	
Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	