


<b>Product:</b>	MediPro Loop	
<b>Trademark:</b>	SAFE®	
<b>REF:</b>	2108	
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG	

<b>Classification:</b>	Class I, based on Medical Device Regulation (EU) 2017/745 MDR, Annex VIII; UMDNS Code: 12-458.
<b>Performance class:</b>	Type IIR (liquid resistant), tested according to EN 14683:2019 for medical face masks
<b>Materials:</b>	<ul style="list-style-type: none"> <li>Outside: liquid-repellent polypropylene fleece</li> <li>Inner side: especially smooth and skin-friendly polypropylene fleece</li> <li>Filter medium: microfiber made of polypropylene</li> <li>Ear loops: soft textile rubber band, free of natural latex</li> <li>Noseband: metal wire, coated with synthetic material</li> </ul>
<b>Product description:</b>	<p>Latex-free, hypoallergenic, high-quality, anatomically shaped face mask with easily modelable nosepiece.</p> <p>The newly developed skin-friendly inner fleece does not lint and, combined with particularly fluffy ear loops, guarantees a high level of wearing comfort.</p> <p>The filter medium consists of special microfiber fleece, which provides an efficient barrier against bacteria.</p> <p>The entire mask is ultrasonically welded using the latest technology.</p> <p>Complies with the HACCP concept.</p>

Product performance:	Product performance according to EN 14683:2019+AC:2019			
	Testing	Performance class	Requirement	Test result
	Bacterial filtration efficiency (BFE), (%)	Type IIR	≥ 98	99,9
	Pressure difference (Pa/cm²)		< 60	24
	Spray resistance pressure (kPa)		≥ 16,0	16,0
	Microbial cleanliness (CFU/g)		≤ 30	≤ 3
Colours:	Green, Blue			

<b>Dermatological compatibility:</b>	No contraindications.	
<b>Packaging and variants:</b>	<b>Packaging</b>	<b>Quantity</b>
	Dispenser box	50 pcs.
	Shipping carton	2.000 pcs.
<b>GTIN:</b>	2108/G	4049825006428
	2108/B	4049825006442
<b>Use:</b>	<p>All materials used are free of irritants, odorless and suitable for medical purposes.</p> <p>Conform to the biological compatibility ISO 10993.</p> <p>The mask is worn by healthcare professionals around the patient to protect the patient.</p>	
<b>Storage:</b>	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 5 years (see labeling on the product).</p>	
<b>Environmental sustainability and disposal:</b>	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>	
<b>Indication:</b>	 <p>Non sterile.</p>	