


Product:	MediCap AstroEco		
Trademark:	SAFE®		
REF:	121G		
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG		
Classification:	Class I of Directive EU 2017/745 on medical devices ; UMDNS code: 15-223.		
Materials:	<ul style="list-style-type: none"><li>▪ Hood: pure non-woven polypropylene, antistatic treated</li><li>▪ Elastic: latex-free synthetic elastomer</li><li>▪ Sewing thread: 100 % polyester</li></ul>		
Product description:	<p>The hood in astronaut shape consists of a head top and a round part.</p> <p>The two parts are joined with an overlock seam.</p> <p>The all-round part can also be used for tying.</p>		
Colour:	Green		
Packaging and Variants:	Packaging	Quantity	
	Bag	100 pcs.	
	Shipping carton	1.000 pcs.	
GTIN:	121G	4049825006404	
Storage:	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 10 years (see labelling on the product).</p>		
Environmental sustainability and disposal:	<p>In the case of contaminated products, the nature and extent of the contamination determines the disposal, following the applicable laws and regulations of the relevant country.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>		