

| | | |
|----------------------|-------------------------------------|---|
| Product: | ShowerCap |  |
| Trademark: | ecoComfort® | |
| REF: | 190 | |
| Manufacturer: | DACH Schutzbekleidung GmbH & Co. KG | |

| | | |
|---|--|-----------------|
| Classification: | CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. | |
| Materials: | <ul style="list-style-type: none"> ▪ Hood: pure polyethylene ▪ Elastic: latex-free synthetic elastomer ▪ Sewing thread: 100 % polyester | |
| Product description: | <p>Waterproof disposable cap in beret shape.</p> <p>The outer finish of the hood is an elastic band bordered with an overlock seam.</p> | |
| Dermatological compatibility: | No intolerance reactions. | |
| Packaging and variants: | Packaging | Quantity |
| | Bag | 100 pcs. |
| | Shipping carton | 1.000 pcs. |
| GTIN: | 190 | 4049825000624 |
| Storage: | <p>Store dry in original packaging without direct sunlight. (see packaging)</p> <p>If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)</p> | |
| Environmental sustainability and disposal: | <p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p> | |