


Product:	EcoCoat	
Trademark:	ecoComfort®	
REF:	410	
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	

Classification:	CAT I according to the Personal Protective Equipment (PPE) Regulation (EU) 2016/425. Class I of Regulation (EU) 2017/745 on medical devices; UMDNS code: 11-897.	
Materials:	<ul style="list-style-type: none"> ▪ Nonwoven fabric: pure polypropylene ▪ Elastic: latex-free, synthetic elastomer ▪ Sewing thread: 100 % polyester 	
Product Description:	EcoCoat is made of soft non-woven fabric, with an opening on the back and two binding belts. The outer finish of the sleeves is elastic edged with an overlock seam. Complies with the HACCP concept.	
Dermatological compatibility:	No intolerance reactions.	
Length:	Size M:115 cm, size L: 136 cm	
Colours:	Green, Blue, White	
Packaging and variants:	Packaging	Quantity
	Bag	10 pcs.
	Shipping carton	100 pcs.
GTIN:	410/B/M	4049825001591
	410/B/L	4049825003885
	410/G/M	4049825001607



	410/G/L	4049825003892
	410/W/M	4049825001614
	410/W/L	4049825003908
Storage:	Store dry in original packaging without direct sunlight. (see packaging) If the storage conditions are observed, the product has a storage life of 10 years. (see labeling on the product)	
Environmental sustainability and disposal:	In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed. A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.	