





879 | GentleSafe® EXTRA | Product Data Sheet

Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG	
Trademark:	SAFE®	
Product:	GentleSafe® EXTRA	
REF:	879/W/S 879/W/M 879/W/L 879/W/XL	

Classification:	CAT III according to Regulation (EU) 2016/425 for personal protective equipment (PPE). Class I according to Regulation (EU) 2017/745 for medical devices.										
Materials:	Carboxylated nitrile butadiene rubber										
Product description:	<p>The nitrile glove is powder-free, latex-free and impervious to liquids.</p> <p>It offers excellent dexterity with high comfort and maximum dexterity (level 5/5). High elasticity with low effort. No hand fatigue.</p> <p>Ambidextrous wear, rolled edge, non-sterile, excellent tactile sensation.</p> <p>Textured surface on the fingertips. Excellent grip, even when handling wet or moist objects.</p> <p>Excellent performance against chemicals, especially cytostatics and solvents, acids and bases, oxidizing agents, viruses, bacteria and fungi.</p> <p>EXTRA longer glove for even more safety against chemicals and EXTRA wearing comfort.</p>										
Pictograms and standards:	<table border="1"> <tr> <td>CAT 3</td> <td>ASTM 1671</td> <td>EN 21420</td> <td>EN 455-2</td> <td>EN ISO 374-5:2016</td> </tr> <tr> <td>AQL1,5</td> <td>ASTM D7160</td> <td>EN 455-1</td> <td>EN 455-3</td> <td>EN ISO 374-1:2016 TYPE C</td> </tr> </table>	CAT 3	ASTM 1671	EN 21420	EN 455-2	EN ISO 374-5:2016	AQL1,5	ASTM D7160	EN 455-1	EN 455-3	EN ISO 374-1:2016 TYPE C
CAT 3	ASTM 1671	EN 21420	EN 455-2	EN ISO 374-5:2016							
AQL1,5	ASTM D7160	EN 455-1	EN 455-3	EN ISO 374-1:2016 TYPE C							

879 | GentleSafe® EXTRA | Product Data Sheet

Other features:	 <p>No intolerance reactions: This product does not contain natural latex.</p>  <p>Virus-tight according to ASTM 1671: Determination of the resistance of materials to the penetration of viral pathogens.</p>  <p>Non sterile.</p>
Length:	min. 300 mm
Width:	S: 85 +/- 5 mm; M: 96 +/- 5 mm; L: 107 +/- 5 mm; XL: 118 +/- 5 mm
Color:	White

Packaging and variants:	Quantity	Packaging
	100 pcs.	Dispenser box
	1000 pcs.	Shipping carton
GTIN:	879/W/S	4049825006688
	879/W/M	4049825006695
	879/W/L	4049825006701
	879/W/XL	4049825006718
Limitation:	<p>The gloves are designed for single use.</p> <p>The duration of protection in working use may differ from the breakthrough time determined according to EN 374, as it depends on the working conditions.</p>	
Storage:	<p>Store in a dry place without direct sunlight in the original packaging (see packaging).</p> <p>If the storage conditions are observed, the product has a storage life of 3 years (see labeling on the product).</p>	

879 | GentleSafe® EXTRA | Product Data Sheet

Environmental sustainability and disposal:	<p>In the case of contaminated products, the type and extent of contamination determines the disposal, and the applicable laws and regulations of the relevant country must be followed.</p> <p>A non-contaminated product can be thermally recycled or disposed of in landfills without releasing toxic substances.</p>
---	--