

Product data sheet

MediCap AstroEco



REF:	122G
Description:	MediCap AstroEco
Trademark:	SAFE®
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG

Categorization:	Class I Medical Device Regulation (EU) 2017/425; UMDNS Code: 15-223.		
Product description:	The bonnet is astronaut shaped and consists of a headpiece and an all round covering piece. Both pieces are connected to each other by an overlock-seam. The all round covering piece can also be tied up. On the forehead covering part a multi layered, absorbent sweatband is added for high comfort and hygiene.		
Material:	Cover:	pure polypropylene nonwoven fabric, antistatic treated	
	Elastic:	latex free, synthetic elastomer	
	Thread:	100% polyester	
Dermatological tolerance:	No intolerance reaction.		
Standard size:	One size.		
Standard color:	Green.		
Packaging unit:	Amounts	Units	GTIN
	100 pcs.	Poly bag	4049825000464
	1000 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 10 years		
Disposal:	This product can be thermally recycled or disposed in a controlled landfill without releasing toxic substances. The disposal of contaminated products must be determined by its contamination. Please respect according laws and regulations.		