

Product data sheet

MediClassic *Excel*




Product:	Surgical Mask
Trademark:	SAFE®
REF:	2101
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG

Categorization:	Class I Medical Device Regulation (EU) 2017/745; UMDNS Code: 12-458. Type II R according to European standard EN 14683:2019 for surgical masks. Note: Type II R is a liquid-resistant type.				
Product description:	Latex free, ergonomically formed face mask with easy adjustable nose piece. The outer layer and the inner layer consist of physiological harmless soft non-woven fabric, the filter media consists of special melt-blown micro fibers, joined by ultrasonic welding.				
Material:	Outer layer:	Fluid repellent polypropylene non-woven fabric			
	Inner layer:	Soft polypropylene non-woven fabric, skin friendly, allergy free			
	Filter media:	Melt-blown micro fiber			
	Strip:	Polypropylene non-woven fabric, skin friendly, allergy free			
	Nose piece:	Polyethylene coated iron wire			
Dermatological tolerance:	All used materials are irritant free and odorless, qualified for medical use.				
Product options:	REF	Color	Units	Packaging	GTIN No.
	2101B	Blue	50 pcs.	Box	
	2101B	Blue	1000 pcs.	Shipping carton	4049825002338
	2101G	Green	50 pcs.	Box	4049825000150
2101G	Green	1000 pcs.	Shipping carton		
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 5 years				



DACH Schutzbekleidung
GmbH & Co. KG

Disposal:	<p>This product can be thermally recycled or disposed in a controlled landfill without releasing toxic substances.</p> <p>The disposal of contaminated products must be determined by its contamination. Please respect according laws and regulations.</p>
Indication:	<p>Non sterile. </p>