

Product data sheet


MediComfort Excel



REF:	221G
Item:	MediComfort Excel
Trademark:	SAFE®
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG

Categorization:	Class I Medical Device Regulation (EU) 2017/745; UMDNS Code: 12-458. Type II R (liquid resistant) according to European standard EN 14683:2019 for surgical masks.		
Product description:	Ergonomically formed face mask with easily adjustable nose piece. The outer layer and the inner layer consist of physiological harmless soft non-woven fabric, the filter media consists of special melt-blown micro fibers, joined by ultrasonic welding.		
Material:	Outer layer: Fluid repellent Polypropylene non-woven fabric Inner layer: Soft Polypropylene non-woven, skin friendly, allergy free Filter media: melt-blown micro fiber Strip: Polypropylene non-woven fabric, skin friendly, allergy free Nose piece: Polyethylene coated iron wire		
Dermatological tolerance:	No intolerance reaction.		
Application:	All used materials are irritant free and odorless, qualified for medical use. This product is suited for permanent usage and provides a brilliant protection.		
Standard color:	Green.		
Packaging unit:	Amounts	Packaging	GTIN 4049825000198
	50 pcs.	Box	
	1000 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life:5 years		



Disposal:	<p>This product can be thermally recycled or disposed in a controlled landfill without releasing toxic substances.</p> <p>The disposal of contaminated products must be determined by its contamination. Please respect according laws and regulations.</p>
Indication:	<p>Non sterile. </p>