

Product data sheet

Stretcher sheet

Product:	Stretcher sheet
Trademark:	ecoComfort®
REF:	3007
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG



Product description:	The stretcher sheet is made of a combination of absorbing cellulose fabric and liquid impermeable foil material. The 17-fold yarn reinforcement guarantees a high tensile strength.		
Material:	Cellulose, polyethylene, nylon thread		
Size / Length:	100 x 210 cm		
Packaging:	Units	Packaging	GTIN No.
	25 pcs.	Polybag	4049825002963
	100 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 10 years		
Disposal:	When used as intended, the product can be disposed of in household waste and landfills and thermally recycled without releasing toxic substances. National regulations and laws must be followed.		