

Product data sheet

# Surgical Gown Easy



|                      |                                     |
|----------------------|-------------------------------------|
| <b>Description:</b>  | Surgical gown                       |
| <b>Trademark:</b>    | SAFE®                               |
| <b>REF:</b>          | 443                                 |
| <b>Manufacturer:</b> | DACH Schutzbekleidung GmbH & Co. KG |

|                                |  |             |              |                  |                 |
|--------------------------------|--|-------------|--------------|------------------|-----------------|
| <b>Categorization:</b>         | Cat I Personal Protective Equipment Regulation (EU) 2016/425.<br>Regulation (EU) 2017/745 on medical devices: Class I  |             |              |                  |                 |
| <b>Essential Requirements:</b> | DIN EN 13795:2011+A1:2013<br>Standard Performance  |             |              |                  |                 |
| <b>Product description:</b>    | The Surgical Gown Easy is made from breathable, abrasion-resistant, and lint-free SMS-material, which consists of three protective layers. With this barrier this surgical gown offers optimal protection against fluids and infective agents, without lowering its wearing comfort or breathability. This gown has its opening on the backside. A Velcro fastener in the neck area and two tie belts in the waist area enable fast dressing and undressing. All seams are ultrasonic welded and thus resistant to fluid and infective agents from penetration. The knitted elastic cuffs ensure a secure fit of gloves. |             |              |                  |                 |
| <b>Material:</b>               | Gown: SMS<br>Outside: fluid repellent polypropylene nonwoven<br>Inside: skin friendly and hypoallergenic polypropylene nonwoven<br>Middle Layer: polypropylene microfiber<br>Knitted cuffs: 100% Polyester<br>Free of PVC (polyvinyl chloride), DEHP(2-ethylhexylphthalat), Latex (natural rubber)<br><b>This product is not sterile.</b>  |             |              |                  |                 |
| <b>Measurement:</b>            | Length: 43 cm, width: 20 cm.   |             |              |                  |                 |
| <b>Product options:</b>        | <b>REF</b>   | <b>Size</b> | <b>Units</b> | <b>Packaging</b> | <b>GTIN No.</b> |
|                                | 443/G/L  | L           | 5 pcs.       | Bag              | 4049825003489   |
|                                | 443/G/L  | L           | 50 pcs.      | Carton           | 4049825003489   |



**DACH Schutzbekleidung**  
GmbH & Co. KG

|                  |   |
|------------------|---|
| <b>Storage:</b>  | Store in a dry environment out of direct sunlight in original packaging, see also packaging.<br>Shelf life: 10 years  |
| <b>Disposal:</b> | When used as intended, the product can be disposed of in household waste and landfills and thermally recycled without releasing toxic substances. National regulations and laws must be followed. |