

Product data sheet

PE Shoe Covers XL

Product:	PE Shoe Covers XL
Trademark:	SAFE®
REF:	952B
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG



Categorization:	Cat I Personal Protective Equipment Regulation (EU) 2016/425.		
Product description:	This overshoe consists of a recyclable synthetic sheet that is liquid proof. The opening is attached with a fused elastic strap.		
Material:	Pure polyethylene Elastic strap: natural rubber		
Size:	Length: 42 cm, width: 25 cm.		
Color:	Blue		
Packaging unit:	Units	Packaging	GTIN No.
	40 pcs.	Poly bag, 5 pieces convoluted	4049825003960
	800 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 10 years		
Disposal:	This product can be thermally recycled or disposed in a controlled landfill without releasing toxic substances. The disposal of contaminated products must be determined by its contamination. Please respect according laws and regulations.		