

Product data sheet

# IsoCoat Zytex® 360



<b>Product:</b>	IsoCoat Zytex® 360
<b>Trademark:</b>	SAFE®
<b>REF:</b>	425Y
<b>Manufacturer:</b>	DACH Schutzbekleidung GmbH & Co. KG

<b>Categorization:</b>	Cat I Personal Protective Equipment Regulation (EU) 2016/425		
<b>Product description:</b>	The isolation gown consists of the material Zytex®, which has been tested according to EN 14126 against the penetration of blood and viruses and according to EN ISO 6529 against the permeation of cytostatics. Zytex® was specially developed for the handling of cytostatics and is highly impermeable to cytostatics. Opening on the back, two tie belts and knitted cuffs.		
<b>Material:</b>	Gown: Zytex® Knitted cuff: 100% polyester Sewing thread: 100% polyester		
<b>Performance:</b>	Class 6: ISO 16603, ASTM F 1670 (resistance to penetration by synthetic blood) Class 6: ISO 16604, ASTM 1671 (resistance to penetration by bacteriophage virus) Class 6: EN 14126, EN ISO 22610 (resistance to bacterial penetration - WET Penetration) Class 3: EN ISO 22612 (resistance to penetration microbial penetration - DRY Penetration) Class 3: ISO 22611 (resistance to penetration by liquid aerosol) Cytostatica resistant: ISO 6529 (resistance to permeation by liquid chemicals)		
<b>Dermatological tolerance:</b>	No intolerance reaction.		
<b>Length:</b>	130 cm		
<b>Color:</b>	Yellow		
<b>Packaging unit:</b>	<b>Units</b>	<b>Packaging</b>	<b>GTIN No.</b>  4049825002444
	10 pcs.	Poly bag	
	50 pcs.	Shipping carton	



**DACH Schutzbekleidung**  
GmbH & Co. KG

<b>Storage:</b>	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 10 years
<b>Disposal:</b>	When used as intended, the product can be disposed of in household waste and landfills and thermally recycled without releasing toxic substances. National regulations and laws must be followed.