

Product data sheet

Inspection Gown



Product:	Inspection Gown
Trademark:	SAFE®
REF:	455
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG

Categorization:	Cat I Personal Protective Equipment Regulation (EU) 2016/425.					
Product description:	The inspection gown consists of soft non-woven fabric, with an open front that can be closed by a velcro fastener. With latex free elastic sleeve straps and overlock stitched seams Complies with the concept of HACCP.					
Material:	Non-woven fabric: pure polypropylene Elastic strap: latex free, synthetic elastomer Thread: 100% polyester					
Dermatological tolerance:	No intolerance reaction.					
Size / length:	Standard size (L): Length 115cm, in the colors blue, red, white, yellow Size S – XL: Length 100 – 125 cm in blue					
Product variants:	REF	Color	Size	Units	Packaging	GTIN No.
	455B/S	Blue	S	10 pcs.	Poly bag	4049825000761
	455B/S	Blue	S	100 pcs.	Shipping carton	
	455B/M	Blue	M	10 pcs.	Poly bag	4049825000754
	455B/M	Blue	M	100 pcs.	Shipping carton	
	455B/L	Blue	L	10 pcs.	Poly bag	4049825000747
	455B/L	Blue	L	100 pcs.	Shipping carton	
455B/XL	Blue	XL	10 pcs.	Poly bag	4049825000778	



	455B/XL	Blue	XL	100 pcs.	Shipping carton	4049825000778
	455G/L	Green	L	10 pcs.	Poly bag	4049825000785
	455G/L	Green	L	100 pcs.	Shipping carton	
	455R/L	Red	L	10 pcs.	Poly bag	4049825000792
	455R/L	Red	L	100 pcs.	Shipping carton	
	455W/L	White	L	10 pcs.	Poly bag	4049825000808
	455W/L	White	L	100 pcs.	Shipping carton	
	455Y/L	Yellow	L	10 pcs.	Poly bag	4049825000815
	455Y/L	Yellow	L	100 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 10 years					
Disposal:	When used as intended, the product can be disposed of in household waste and landfills and thermally recycled without releasing toxic substances. National regulations and laws must be followed.					